


*Hazardous according to criteria of Australian Safety and Compensation Council.*

## 1 Identification

- **Product identifier**
- **Trade name:** IVISbrite® D-Luciferin Potassium Salt Bioluminescent Substrate
- **Product number:** 122799, 122799-10, 122799-5, 122799-100MG
- **CAS Number:**  
115144-35-9
- **Relevant identified uses of the substance or mixture and uses advised against**
- **Product category** PC21 Laboratory chemicals
- **Application of the substance / the mixture** Laboratory chemicals
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
Supplier/Local:  
-----  
Revvity Pty. Ltd  
Building C, Level 2, Tenancy A,  
211 Wellington Road  
Mulgrave 3170, VIC Australia
- **Emergency telephone number:**  
CHEMTREC (within US) 800-424-9300  
CHEMTREC (from outside US) 1(703)-527-3887
- **Local Emergency Number CHEMTREC - (Within Australia) +(61)-290372994**  
Revvity, Inc  
549 Albany Street  
Boston, MA 02118
- **Further information obtainable from:**  
US Technical Support  
800-762-4000

## 2 Hazard(s) Identification

- **Classification of the substance or mixture**  
*Serious eye damage/irritation – Category 2A H319 Causes serious eye irritation.*
  - **Additional information:** *For the wording of the relevant risk phrases refer to section 16.*
  - 
  - **Label elements**
  - **GHS label elements**  
*The substance is classified and labelled according to the Globally Harmonised System (GHS).*
  - **Hazard pictograms**
- 

GHS07
- **Signal word** *Warning*
  - **Hazard statements**  
*Causes serious eye irritation.*
  - **Precautionary statements**  
*Wash thoroughly after handling.*

(Contd. on page 2)

# Safety Data Sheet

## according to WHS Regulations

Printing date 11.03.2024

Revision: 31.07.2023

**Trade name: IVISbrite® D-Luciferin Potassium Salt Bioluminescent Substrate**

(Contd. of page 1)

*Wear eye protection / face protection.*

*IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.*

*Continue rinsing.*

*If eye irritation persists: Get medical advice/attention.*

- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### 3 Composition and Information on Ingredients

- **Chemical characterisation: Substances**
- **CAS No. Description**  
115144-35-9 D-Luciferin, Potassium Salt

### 4 First Aid Measures

- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Immediately rinse with water.
- **After eye contact:**  
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:** If symptoms persist consult doctor.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### 5 Fire Fighting Measures

- **Suitable extinguishing agents:** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Protective equipment:** Wear self-contained respiratory protective device.

### 6 Accidental Release Measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:** Pick up mechanically.
- **Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and Storage

- **Handling:**
- **Precautions for safe handling**  
Store in cool, dry place in tightly closed receptacles.  
No special precautions are necessary if used correctly.

(Contd. on page 3)

# Safety Data Sheet

## according to WHS Regulations

Printing date 11.03.2024

Revision: 31.07.2023

**Trade name: IVISbrite® D-Luciferin Potassium Salt Bioluminescent Substrate**

(Contd. of page 2)

- **Information about fire - and explosion protection:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and containers:** No special requirements.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep container tightly sealed.
- **Storage class:** 11
- **Specific end use(s)** No further relevant information available.

### 8 Exposure controls and personal protection

- **Additional information about design of technical facilities:** No further data; see section 7.
- **Ingredients with limit values that require monitoring at the workplace:** Not required.
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Keep away from foodstuffs, beverages and feed.  
Immediately remove all soiled and contaminated clothing  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes.  
Avoid contact with the eyes and skin.
- **Respiratory protection:** Suitable respiratory protective device recommended.
- **Protection of hands:**



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**  
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:**



Tightly sealed goggles

### 9 Physical and Chemical Properties

- **General Information**
- **Appearance:**
- **Form:** Powder
- **Colour:** Whitish
- **Odour:** Odourless
- **Odour threshold:** Not determined.
- **pH-value:** Not applicable.
- **Change in condition**
- **Melting point/freezing point:** Undetermined.

(Contd. on page 4)

-AU

# Safety Data Sheet

## according to WHS Regulations

Printing date 11.03.2024

Revision: 31.07.2023

**Trade name: IVISbrite® D-Luciferin Potassium Salt Bioluminescent Substrate**

(Contd. of page 3)

- **Initial boiling point and boiling range:** Undetermined.
- **Flash point:** Not applicable.
- **Flammability (solid, gas):** Product is not flammable.
- **Decomposition temperature:** Not determined.
- **Ignition temperature:** Not determined.
- **Explosive properties:** Product does not present an explosion hazard.
- **Explosion limits:**
- **Lower:** Not determined.
- **Upper:** Not determined.
- **Vapour pressure:** Not applicable.
- **Density:** Not determined.
- **Relative density** Not determined.
- **Vapour density** Not applicable.
- **Evaporation rate** Not applicable.
- **Solubility in / Miscibility with**
- **water:** Insoluble.
- **Partition coefficient: n-octanol/water:** Not determined.
- **Viscosity:**
- **Dynamic:** Not applicable.
- **Kinematic:** Not applicable.
- **Other information** No further relevant information available.

### 10 Stability and Reactivity

- **Reactivity** No further relevant information available.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

### 11 Toxicological Information

- **Information on toxicological effects**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Causes serious eye irritation.
- **Respiratory or skin sensitisation** Based on available data, the classification criteria are not met.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.
- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

### 12 Ecological Information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.

(Contd. on page 5)

## Safety Data Sheet according to WHS Regulations

Printing date 11.03.2024

Revision: 31.07.2023

**Trade name: IVISbrite® D-Luciferin Potassium Salt Bioluminescent Substrate**

(Contd. of page 4)

- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

### 13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**  
Must not be disposed together with household garbage. Do not allow product to reach sewage system.  
Must be specially treated adhering to official regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

### 14 Transport information

· <b>UN-Number</b> · <b>ADG, IMDG, IATA</b>	Void
· <b>UN proper shipping name</b> · <b>ADG, IMDG, IATA</b>	Void
· <b>Transport hazard class(es)</b> · <b>ADG, ADN, IMDG, IATA</b> · <b>Class</b>	Void
· <b>Packing group</b> · <b>ADG, IMDG, IATA</b>	Void
· <b>Environmental hazards:</b>	Not applicable.
· <b>Special precautions for user</b>	Not applicable.
· <b>Transport in bulk according to Annex II of Marpol and the IBC Code</b>	Not applicable.
· <b>UN "Model Regulation":</b>	Void

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

· <b>Australian Inventory of Industrial Chemicals</b> Substance is not listed.
· <b>Standard for the Uniform Scheduling of Medicines and Poisons</b> Substance is not listed.
· <b>Australia: Priority Existing Chemicals</b> Substance is not listed.

(Contd. on page 6)

AU

# Safety Data Sheet

## according to WHS Regulations

Printing date 11.03.2024

Revision: 31.07.2023

**Trade name: IVISbrite® D-Luciferin Potassium Salt Bioluminescent Substrate**

(Contd. of page 5)

- **GHS label elements**

The substance is classified and labelled according to the Globally Harmonised System (GHS).

- **Hazard pictograms**



GHS07

- **Signal word** Warning

- **Hazard statements**

Causes serious eye irritation.

- **Precautionary statements**

Wash thoroughly after handling.

Wear eye protection / face protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If eye irritation persists: Get medical advice/attention.

- **Directive 2012/18/EU**

· **Named dangerous substances - ANNEX I** Substance is not listed.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

## 16 Other information

The information provided in this safety data sheet is based on our current knowledge, and is believed to be correct at the date of publication. However, no representation is made concerning its accuracy and completeness. It is intended as guidance only, and is not to be regarded as a warranty or specification of quality. All materials may present unknown hazards and should be used with caution. Although certain hazards are described, we cannot guarantee that these are the only hazards that exist. Revvity, Inc. cannot be held liable for any damage resulting from handling or contact with the product.

- **Contact:**

- **Abbreviations and acronyms:**

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Serious eye damage/irritation – Category 2A: Serious eye damage/eye irritation – Category 2A

AU